

GLP-2 (1-34) (Human) - Biotin Labeled Purified IgG

Catalog#	Standard size	Price	Order
B-G-028-14	100 µl	\$450.00	<input type="button" value="ORDER"/>

Product description

More information

References

[GLP-1 & GLP-2](#)

[Normal Plasma Peptide/Protein Levels](#)

[Neuron Function Modulation](#)

Catalog #	B-G-028-14
Standard Size	100 µl
Sequence	His - Ala - Asp - Gly - Ser - Phe - Ser - Asp - Glu - Met - Asn - Thr - Ile - Leu - Asp - Asn - Leu - Ala - Ala - Arg - Asp - Phe - Ile - Asn - Trp - Leu - Ile - Gln - Thr - Lys - Ile - Thr - Asp - Arg
Species	Human
Storage Condition	Store at 4°C
Content	This vial contains 100µl of rabbit anti-Glucagon-Like Peptide 2 (GLP-2) (Human) IgG biotinylated in PBS (137mM NaCl, 2.7mM KCl, 4.3mM Na ₂ HPO ₄ , 1.4mM KH ₂ PO ₄ , pH 7.4), which was purified form whole serum with Portein A Column.

Thursday 20 December, 2012

11983127 requests since Wednesday 05 April, 2006

Copyright © 2012 PHOENIX PHARMACEUTICALS, INC.