

Glicentin-related Polypeptide (21-50) (Human) - Biotin Labeled Purified IgG

Catalog#	Standard size	Price	Order
B-G-028-17	100 µl	\$500.00	<input type="button" value="ORDER"/>

Product description

References

Catalog #	B-G-028-17
Standard Size	100 µl
Sequence	Arg - Ser - Leu - Gln - Asp - Thr - Glu - Glu - Lys - Ser - Arg - Ser - Phe - Ser - Ala - Ser - Gln - Ala - Asp - Pro - Leu - Ser - Asp - Pro - Asp - Gln - Met - Asn - Glu - Asp
Species	Human
Storage Condition	The antibody can be stored at 4 °C for one or two weeks without detectable loss of activity. Antibody can also be aliquotted and stored frozen at -20°C to -70 °C. in a manual defrost freezer for six months. Avoid repeat freeze-thaw of this product for best results.
Content	This vial contains 100µl of rabbit IgG biotinylated in PBS (137mM NaCl, 2.7mM KCl, 4.3mM Na ₂ HPO ₄ , 1.4mM KH ₂ PO ₄ , pH 7.4), which was purified form whole serum with Portein A Column.

Thursday 20 December, 2012

11983131 requests since Wednesday 05 April, 2006

Copyright © 2012 PHOENIX PHARMACEUTICALS, INC.