

Exendin-4 (Heloderma suspectum) - Biotin Labeled Purified IgG

Catalog#	Standard size	Price	Order
B-G-070-94	100 µl	\$500.00	<input type="button" value="ORDER"/>

Product description

More information

References

[Exendin-4, Research Abstracts, Kits Standard Curves, Sequences, etc.](#)

[Neuron Function Modulation](#)

Catalog #	B-G-070-94
Standard Size	100 µl
Sequence	His - Gly - Glu - Gly - Thr - Phe - Thr - Ser - Asp - Leu - Ser - Lys - Gln - Met - Glu - Glu - Glu - Ala - Val - Arg - Leu - Phe - Ile - Glu - Trp - Leu - Lys - Asn - Gly - Gly - Pro - Ser - Ser - Gly - Ala - Pro - Pro - Pro - Ser - NH ₂
Storage Condition	Store at 4°C
Content	This vial contains biotin labeled rabbit IgG in 100µl PBS (137mM NaCl, 2.7mM KCl, 4.3mM Na ₂ HPO ₄ , 1.4mM KH ₂ PO ₄ , pH 7.4) in the lyophilized form.

Thursday 20 December, 2012

11983231 requests since Wednesday 05 April, 2006

Copyright © 2012 PHOENIX PHARMACEUTICALS, INC.