

# Oxyntomodulin (Human, Rat, Mouse) - EIA Kit

Catalog#	Standard size	Price	Order
EK-028-22	1 kit	\$420.00	<input type="button" value="ORDER"/>

## Product description

## More information

Oxyntomodulin

Oxyntomodulin (Human, Rat, Mouse) - EIA Kit

<b>Protocol</b>	Assay Protocol (PDF) *
<b>Sample Preparation</b>	Blood collection Plasma extraction Tissue extraction CSF extraction
<b>Catalog #</b>	EK-028-22
<b>Standard Size</b>	96 Wells
<b>Sensitivity</b>	Minimum Detectable Concentration = 0.25 ng/ml
<b>Range</b>	0-100 ng/ml
<b>Linear Range</b>	0.25-3.9 ng/ml

Cross Reactivity	Peptides	Cross-reactivity
	Oxyntomodulin (Human, Rat, Mouse)	100
	Oxyntomodulin (Porcine)	100
	Oxyntomodulin (30-37) (Porcine)	3
	GLP-1 (7-37) (Human, Rat, Mouse)	0
	GLP-1 (7-36)-amide (Human, Rat, Mouse)	0
	GLP-2 (Human)	0
	Glucagon (Human)	0
	Insulin (Human)	0

**Sequence** His - Ser - Gln - Gly - Thr - Phe - Thr - Ser - Asp - Tyr - Ser - Lys - Tyr - Leu - Asp - Ser - Arg - Arg - Ala - Gln - Asp - Phe - Val - Gln - Trp - Leu - Met - Asn - Thr - Lys - Arg - Asn - Arg - Asn - Asn - Ile - Ala

Standard Curve Image



