

OBAP (Rat, Mouse) - Cy3 Labeled

| Catalog# | Standard size | Price | Order |
|------------|---------------|----------|--------------------------------------|
| FC3-055-51 | 1 nmol | \$450.00 | <input type="button" value="ORDER"/> |

Product description

More information

Osteoblast Activating Peptide (OBAP)

| | |
|--------------------------|---|
| Catalog # | FC3-055-51 |
| Standard Size | 1 nmol |
| Sequence | Cy3 - (Leu - Asp - Leu - Asn - Leu - Asp - Leu - Ser - Lys - Phe - Arg - Leu - Pro - Gln - Pro - Ser - Ser - Gly - Arg - Glu - Ser - Pro - Arg - His) |
| Species | Rat, Mouse |
| M.W | 3377.4 |
| Solubility | Soluble in 50mM sodium phosphate buffer (pH 8.5) or 10% DMSO in PBS buffer (pH 7.5). |
| Storage Condition | Store at -20°C with desiccant. |
| Contents | The contents of this vial have been accurately determined to contain 1 nmole of NET peptide. |
| Absorption | Excitation maximum: 558 nm Extinction coefficient: 130000 M ⁻¹ cm ⁻¹ |
| Fluorescent_Dye | Cy TM 3 |
| Emission_maximum | 572 nm |

Thursday 20 December, 2012

11983669 requests since Wednesday 05 April, 2006

Copyright © 2012 PHOENIX PHARMACEUTICALS, INC.