

OBAP (Human) - Cy3 Labeled

Catalog#	Standard size	Price	Order
FC3-055-52	1 nmol	\$450.00	<input type="button" value="ORDER"/>

Product description

More information

Osteoblast Activating Peptide (OBAP)

Catalog #	FC3-055-52
Standard Size	1 nmol
Sequence	Cy3 - (Leu - Asp - Leu - Asn - Leu - Glu - Leu - Ser - Lys - Phe - Arg - Met - Pro - Gln - Pro - Ser - Ser - Gly - Arg - Glu - Ser - Pro - Arg - His)
Species	Human
M.W	3409.4
Solubility	Soluble in 50mM sodium phosphate buffer (pH 8.5) or 10% DMSO in PBS buffer (pH 7.5).
Storage Condition	Store at -20°C with desiccant.
Contents	The contents of this vial have been accurately determined to contain 1 nmole of NET peptide.
Absorption	Excitation maximum: 558 nm Extinction coefficient: 130000 M ⁻¹ cm ⁻¹
Fluorescent_Dye	Cy TM 3
Emission_maximum	572 nm