

# Adipo R2 (374-386) / Adiponectin Receptor Type 2 (374-386) (Human) - Cy3 Labeled Purified IgG

Catalog#	Standard size	Price	Order
FC3-G-001-23	100 µl	\$650.00	<input type="button" value="ORDER"/>

## Product description

## More information

## References

[Adiponectin, Adiponectin Receptors & Adiponectin EIA Kits](#)

[Adiponectin Receptor-1 &-2 \(AdipoR1 & AdipoR2\) two recently identified receptors](#)

[Normal Plasma Peptide/Protein Levels](#)

[Predict Alzheimer's Disease up to 10 years in advance](#)

<b>Catalog #</b>	FC3-G-001-23
<b>Standard Size</b>	100 µl
<b>Sequence</b>	Phe - Met - Ile - Gly - Gly - Gly - Cys - Ser - Glu - Glu - Asp - Ala - Leu
<b>Species</b>	Human
<b>Storage Condition</b>	Upon receipt, store product in its original form at -20°C (or below) and keep away from light, especially after reconstitution.
<b>Contents</b>	Each vial contains 100µl of lyophilized Cy3 labeled rabbit anti-Adiponectin Receptor-2 (374-386) [Adipo R2 (374-386)] (Human) IgG.
<b>Absorption</b>	Excitation maximum: 558 nm Extinction coefficient: 150000 M <sup>-1</sup> cm <sup>-1</sup>
<b>Fluorescent_Dye</b>	Cy <sup>TM</sup> 3
<b>Emission_maximum</b>	572 nm

Thursday 20 December, 2012

11983691 requests since Wednesday 05 April, 2006

Copyright © 2012 PHOENIX PHARMACEUTICALS, INC.