

Adipo R1 (357-375) / Adiponectin Receptor Type 1 (357-375) (Human) - Cy3 Labeled Purified IgG

Catalog#	Standard size	Price	Order
FC3-G-001-44	100 µl	\$600.00	<input type="button" value="ORDER"/>

Product description

More information

References

[Adiponectin, Adiponectin Receptors & Adiponectin EIA Kits](#)

[Adiponectin Receptor-1 &-2 \(AdipoR1 & AdipoR2\) two recently identified receptors](#)

[Normal Plasma Peptide/Protein Levels](#)

[Predict Alzheimer's Disease up to 10 years in advance](#)

Catalog #	FC3-G-001-44
Standard Size	100 µl
Sequence	Asn - Leu - Gln - Glu - Phe - Arg - Tyr - Gly - Leu - Glu - Gly - Gly - Cys - Thr - Asp - Asp - Thr - Leu - Leu
Species	Human
Storage Condition	Upon receipt, store product in its original form at -20°C (or below) and keep away from light, especially after reconstitution.
Contents	Each vial contains 100ul of lyophilized Cy3 labeled rabbit anti-Adiponectin Receptor-1 (357-375)-NH2 [Adipo R1 (357-375)- NH2] (Human) IgG.
Absorption	Excitation maximum: 558 nm Extinction coefficient: 150000 M ⁻¹ cm ⁻¹
Fluorescent_Dye	Cy TM 3
Emission_maximum	572 nm

Thursday 20 December, 2012

11983692 requests since Wednesday 05 April, 2006

Copyright © 2012 PHOENIX PHARMACEUTICALS, INC.