

Growth Hormone Secretagogue Receptor type 1B (245-289) (Human) - Cy3 Labeled Purified IgG

Catalog#	Standard size	Price	Order
FC3-G-001-61	100 µl	\$550.00	<input type="button" value="ORDER"/>

Product description

References

Catalog #	FC3-G-001-61
Standard Size	100 µl
Sequence	Gly - Asp - Ala - Val - Val - Gly - Ala - Ser - Leu - Arg - Asp - Gln - Asn - His - Lys - Gln - Thr - Val - Lys - Met - Leu - Gly - Gly - Ser - Gln - Arg - Ala - Leu - Arg - Leu - Ser - Leu - Ala - Gly - Pro - Ile - Leu - Ser - Leu - Cys - Leu - Leu - Pro - Ser - Leu
Species	Human
Storage Condition	Upon receipt, store product in its original form at -20°C (or below) and keep away from light, especially after reconstitution.
Contents	Each vial contains 100µl of lyophilized Cy3 labeled rabbit anti- Growth Hormones Secretagogue Receptor Type 1b (245-289) (Human) IgG.
Absorption	Excitation maximum: 558 nm Extinction coefficient: 150000 M ⁻¹ cm ⁻¹
Fluorescent_Dye	Cy TM 3
Emission_maximum	572 nm

Thursday 20 December, 2012

11983702 requests since Wednesday 05 April, 2006

Copyright © 2012 PHOENIX PHARMACEUTICALS, INC.