

Visfatin (412-431) (Human) - Cy3 Labeled Purified IgG

Catalog#	Standard size	Price	Order
FC3-G-003-82	100 µl	\$650.00	<input type="button" value="ORDER"/>

Product description


More information

References

Visfatin

Normal Plasma Peptide/Protein Levels

Catalog #	FC3-G-003-82
Standard Size	100 µl
Sequence	Asn - Val - Phe - Lys - Asp - Pro - Val - Ala - Asp - Pro - Asn - Lys - Arg - Ser - Lys - Lys - Gly - Arg - Leu - Ser
Species	Human
Storage Condition	Upon receipt, store product in its original form at -20°C (or below) and keep away from light, especially after reconstitution.
Contents	Each vial contains 100µl of Cy3 labeled rabbit anti-Visfatin (412-431) (Human) IgG.
Absorption	Absorbance Max: 550 nm Extinction Max: 150000 M ⁻¹ cm ⁻¹
Fluorescent_Dye	Cy TM 3
Emission_maximum	572 nm
Image	1/200x .1/100x .1/50x



Serial dilution of labeled antibody scanned by fluorescence scanner at 532nm (4100A, Axon Instrument).