

Humanin, Mitochondria-encoded / HN(M) (Human) - Cy3 Labeled Purified IgG

| Catalog# | Standard size | Price | Order |
|--------------|---------------|----------|--------------------------------------|
| FC3-G-018-25 | 100 µl | \$550.00 | <input type="button" value="ORDER"/> |

Product description

More information

Neuron Function Modulation

| | |
|----------------------|--------------|
| Catalog # | FC3-G-018-25 |
| Standard Size | 100 µl |
| Species | Human |

Thursday 20 December, 2012

11983764 requests since Wednesday 05 April, 2006

Copyright © 2012 PHOENIX PHARMACEUTICALS, INC.