

## MGF (Human) - Cy3 Labeled Purified IgG

Catalog#	Standard size	Price	Order
FC3-G-033-35	100 µl	\$650.00	<input type="button" value="ORDER"/>

### Product description

### More information

### References

[MGF \(mechano-growth factor\)](#)

[FGL Peptide and BCL Two Ligands of FGFR](#)

<b>Catalog #</b>	FC3-G-033-35
<b>Standard Size</b>	100 µl
<b>Sequence</b>	Tyr - Gln - Pro - Pro - Ser - Thr - Asn - Lys - Asn - Thr - Lys - Ser - Gln - Arg - Arg - Lys - Gly - Ser - Thr - Phe - Glu - Glu - Arg - Lys
<b>Species</b>	Human
<b>Storage Condition</b>	Upon receipt, store product in its original form at -20°C (or below) and keep away from light, especially after reconstitution.
<b>Contents</b>	Each vial contains 100µl of lyophilized Cy3 labeled rabbit Anti- Mechano Growth Factor (MGF) / C-Terminal Peptide of IGF-1 Ec (Human) IgG.
<b>Absorption</b>	Excitation maximum: 558 nm Extinction coefficient: 150000 M <sup>-1</sup> cm <sup>-1</sup>
<b>Fluorescent_Dye</b>	Cy <sup>TM</sup> 3
<b>Emission_maximum</b>	572 nm

Thursday 20 December, 2012

11983799 requests since Wednesday 05 April, 2006

Copyright © 2012 PHOENIX PHARMACEUTICALS, INC.