

## MSH, Alpha (Human, Rat, Mouse) - Cy3 Labeled Purified IgG

Catalog#	Standard size	Price	Order
FC3-G-043-01	100 µl	\$595.00	<input type="button" value="ORDER"/>

### Product description

### More information

### References

#### Alpha-MSH (Alpha-Melanocyte Stimulating Hormone)

<b>Catalog #</b>	FC3-G-043-01
<b>Standard Size</b>	100 µl
<b>Sequence</b>	Ac - Ser - Tyr - Ser - Met - Glu - His - Phe - Arg - Trp - Gly - Lys - Pro - Val - NH <sub>2</sub>
<b>Species</b>	Human, Rat, Mouse
<b>Storage Condition</b>	Upon receipt, store product in its original form at -20°C (or below) and keep away from light.
<b>Contents</b>	Each vial contains 100µl of Cy3 labeled rabbit anti-alpha-MSH (Human, Rat, Mouse) IgG.
<b>Absorption</b>	Absorbance Max: 550 nm Extinction Max: 150000 M <sup>-1</sup> cm <sup>-1</sup>
<b>Fluorescent_Dye</b>	Cy <sup>TM</sup> 3
<b>Emission_maximum</b>	570 nm
<b>Image</b>	1/400x ./.200x ./.100x



Serial dilution of labeled antibody scanned by fluorescence scanner (4100A, Axon Instrument) using a 532nm laser