

Xenin 25 (Human, Rat, Mouse) - Cy3 Labeled Purified IgG

Catalog#	Standard size	Price	Order
FC3-G-046-74	100 µl	\$475.00	<input type="button" value="ORDER"/>

Product description

More information

Xenin-25

Catalog #	FC3-G-046-74
Standard Size	100 µl
Sequence	Met - Leu - Thr - Lys - Phe - Glu - Thr - Lys - Ser - Ala - Arg - Val - Lys - Gly - Leu - Ser - Phe - His - Pro - Lys - Arg - Pro - Trp - Ile - Leu
Species	Human, Rat, Mouse

Thursday 20 December, 2012

11983816 requests since Wednesday 05 April, 2006

Copyright © 2012 PHOENIX PHARMACEUTICALS, INC.