

Pancreatic Polypeptide (Human) - Cy5 Labeled

Catalog#	Standard size	Price	Order
FC5-054-02	1 nmol	\$400.00	<input type="button" value="ORDER"/>

- Product description
- More information
- References

[Predict Alzheimer's Disease up to 10 years in advance](#)

Catalog #	FC5-054-02
Standard Size	1 nmol
Species	Human
M.W	4859.77
Solubility	Soluble in 50mM sodium phosphate buffer (pH 8.5) or 10% DMSO in PBS buffer (pH 7.5).
Storage Condition	Store at -20°C with desiccant.
Contents	This contents of this vial have been accurately determined to contain 600µg of NET peptide
Absorption	Excitation maximum: 646 nm Extinction coefficient: 250000 M ⁻¹ cm ⁻¹
Fluorescent_Dye	Cy TM 5
Emission_maximum	664 nm