

Apelin-15 (63-75) (Human, Rat, Mouse, Bovine) - Cy5 Labeled

Catalog#	Standard size	Price	Order
FC5-057-27	1 nmol	\$250.00	<input type="button" value="ORDER"/>

- Product description
- More information
- References

[Apelin - Research Abstracts, Antibody Stainings, Kits Standard Curve, Sequences, etc.](#)

[Normal Plasma Peptide/Protein Levels](#)

Catalog #	FC5-057-27
Standard Size	1 nmol
Sequence	Cy5 - (Arg - Arg - Gln - Arg - Pro - Arg - Leu - Ser - His - Lys - Gly - Pro - Met)
Species	Human, Bovine, Rat, Mouse
M.W	2258.87
Solubility	Soluble in 50mM sodium phosphate buffer (pH 8.5) or 10% DMSO in PBS buffer (pH 7.5).
Storage Condition	Store at -20°C with desiccant.
Contents	The contents of this vial have been accurately determined to contain 1 nmole of NET peptide
Absorption	Excitation maximum: 646 nm Extinction coefficient: 250000 M ⁻¹ cm ⁻¹
Fluorescent_Dye	Cy TM 5
Emission_maximum	664 nm