


Growth Hormone Secretagogue Receptor type 1A [Cys0] (330-366) (Human) - Cy5 Labeled Purified IgG

Catalog#	Standard size	Price	Order
FC5-G-001-62	100 µl	\$550.00	<input type="button" value="ORDER"/>

Product description

References

Catalog #	FC5-G-001-62
Standard Size	100 µl
Sequence	Cys - Tyr - Arg - Val - Ala - Val - Phe - Arg - Leu - Leu - Gly - Phe - Glu - Pro - Phe - Ser - Gln - Arg - Lys - Leu - Ser - Thr - Leu - Lys - Asp - Glu - Ser - Ser - Arg - Ala - Trp - Thr - Glu - Ser - Ser - Ile - Asn - Thr
Species	Human
Storage Condition	Upon receipt, store at -20°C (or below) and keep away from light.
Contents	Each vial contains 100µl of Cy5 labeled rabbit anti- [Cys0]-GHS-R type 1a (330-366) (Human) purified IgG in liquid form.
Absorption	Absorbance Max: 649 nm Extinction Max: 250 000 M ⁻¹ cm ⁻¹
Fluorescent_Dye	Cy TM 5
Emission_maximum	670 nm
Image	62.5 125 250 500 1000 µg/ml



Serial peptide dilutions (25nl/spot) detected by 1:100 dilution of labeled antibody scanned by fluorescence scanner at 635nm (4100A, Axon Instrument).