

# Visfatin (400-450) (Human) - Cy5 Labeled Purified IgG

Catalog#	Standard size	Price	Order
FC5-G-003-84	100 µl	\$650.00	<input type="button" value="ORDER"/>

## Product description

## More information

## References

### Visfatin

#### Normal Plasma Peptide/Protein Levels

<b>Catalog #</b>	FC5-G-003-84
<b>Standard Size</b>	100 µl
<b>Sequence</b>	Lys - Cys - Ser - Tyr - Val - Val - Thr - Asn - Gly - Leu - Gly - Ile - Asn - Val - Phe - Lys - Asp - Pro - Val - Ala - Asp - Pro - Asn - Lys - Arg - Ser - Lys - Lys - Gly - Arg - Leu - Ser - Leu - His - Arg - Thr - Pro - Ala - Gly - Asn - Phe - Val - Thr - Leu - Glu - Glu - Gly - Lys - Gly - Asp - Leu
<b>Species</b>	Human
<b>Storage Condition</b>	Upon receipt, store product in its original form at -20°C (or below) and keep away from light, especially after reconstitution.
<b>Contents</b>	Each vial contains 100µl of Cy5 labeled rabbit anti-Visfatin (400-450) (Human) IgG.
<b>Absorption</b>	Absorbance Max: 649 nm Extinction Max: 250 000 M <sup>-1</sup> cm <sup>-1</sup>
<b>Fluorescent_Dye</b>	Cy <sup>TM</sup> 3
<b>Emission_maximum</b>	670 nm
<b>Image</b>	1/200x . /100x . /50x



Serial dilution of labeled antibody scanned by fluorescence scanner at 635nm (4100A, Axon Instrument).