

# TLQP-21 (Rat, Mouse), Fluorescent EIA Kit

Catalog#	Standard size	Price	Order
FEK-003-89	1 kit	\$460.00	<input type="button" value="ORDER"/>

**Product description**    **More information**

TLQP and other VGF-derived Peptides  
 Orphan Bioactive Peptides from VGF and PC2  
 Alzheimer biomarker

<b>Protocol</b>	Fluorescent EIA Kit (PDF)
<b>Sample Preparation</b>	Blood collection Plasma extraction Tissue extraction CSF extraction
<b>Catalog #</b>	FEK-003-89
<b>Standard Size</b>	96 Well Kit
<b>Sensitivity</b>	Minimum Detectable Concentration = 34.8 pg/ml
<b>Range</b>	0~10000 pg/ml
<b>Linear Range</b>	34.8 ~ 550 pg/ml

Cross Reactivity	Peptide name	% Cross-reactivity
	TLQP-21 (Rat, Mouse)	100
	TLQP-21 (Human)	61
	TLQP-21 (Rat)	29
	alpha-ANP (Human, Rat, Mouse, Porcine)	0
	AVP (Human, Rat, Mouse, Porcine)	0
	BNP-32 (Human)	0
	CNP-22 (Human, Rat, Mouse, Porcine)	0
	VIP (Human,Rat,Porcine)	0
	Adrenomedullin(1-52)(Human)	0
	PACAP 38 (Human, Rat, Ovine)	0
	Endothelin-1 (Human, Rat, Mouse, Canine, Bovine, Porcine)	0
	Prepro-VGF (586-615) /AQEE-30 (Human)	0

Standard Curve Image



