

TLQP-21 (Human) - Fluorescent EIA Kit

Catalog#	Standard size	Price	Order
FEK-003-90	1 kit	\$460.00	ORDER

Product description

More information

[TLQP and other VGF-derived Peptides](#)

[Orphan Bioactive Peptides from VGF and PC2](#)

[Alzheimer biomarker](#)

Protocol [Fluorescent EIA Kit \(PDF\)](#)

Sample Preparation [Blood collection](#)
[Plasma extraction](#)
[Tissue extraction](#)
[CSF extraction](#)

Catalog # FEK-003-90

Standard Size 96 well


Sensitivity Minimum Detectable Concentration = 21.4 pg/ml

Range 0-10000 pg/ml

Linear Range 21.4-381 pg/ml

Cross Reactivity	Peptide	%Cross-reactivity
	TLQP-21 (Human)	100
	TLQP-21 (Rat)	30.2
	TLQP-62 (Rat)	0
	alpha-Atrial Natriuretic Polypeptide (alpha-ANP) (Human)	0
	[Arg8]-Vasopressin (AVP) (Human, Rat, Mouse, Ovine)	0
	C-Type Natriuretic Peptide-22 (CNP-22) (Human)	0
	Brain Natriuretic Peptide-32 (BNP-32) (Human)	0
	Prepro-VGF (586-615) / AQEE-30 (Human)	0
	Vasoactive Intestinal Peptide (VIP)	0
	Adrenomedullin (1-52)(Human)	0
	PACAP 38 (Human, Rat, Ovine)	0
	Endothelin-1 (Human)	0

Standard
Curve Image



4-P Fit: $y = (A - D) / (1 + (x/C)^B) + D$