

APL 1beta-28 - Fluorescent EIA Kit

Catalog#	Standard size	Price	Order
FEK-018-42	1 kit	\$460.00	<input type="button" value="ORDER"/>

Product description

More information

[Amyloid-beta Precursor-like Proteins \(APLPs\)](#)

[Peptide Fingerprint of Alzheimer s](#)

[Alzheimer's Disease fingerprint Peptide Library](#)

Protocol [Fluorescent EIA Kit \(PDF\)](#)

Sample Preparation [Blood collection](#)
[Plasma extraction](#)
[Tissue extraction](#)
[CSF extraction](#)

Catalog # FEK-018-42

Standard Size 96 wells

Sensitivity Minimum Detectable Concentration = 51 pg/ml

Range 0-10000 pg/ml

Linear Range 51 - 1125 pg/ml

Cross Reactivity	Peptide	%Cross-reactivity
	APL 1 β -28	100
	APL 1 β -27	72
	APL 1 β -28 (1-14)	0
	Amyloid b-Protein (1-40) (Human)	0
	Amyloid b-Protein (1-42)	0
	Amyloid b-Protein (1-43)	0

Sequence Asp - Glu - Leu - Ala - Pro - Ala - Gly - Thr - Gly - Val - Ser - Arg - Glu - Ala - Val - Ser - Gly - Leu - Leu - Ile - Met - Gly - Ala - Gly - Gly - Gly - Ser - Leu

Standard
Curve Image

