

Growth Hormone Secretagogue Receptor type 1B (245-289) (Human) - Rhodamine Labeled Purified IgG

Catalog#	Standard size	Price	Order
FR-G-001-61	100 µl	\$450.00	<input type="button" value="ORDER"/>

Product description

References

Catalog #	FR-G-001-61
Standard Size	100 µl
Sequence	Gly - Asp - Ala - Val - Val - Gly - Ala - Ser - Leu - Arg - Asp - Gln - Asn - His - Lys - Gln - Thr - Val - Lys - Met - Leu - Gly - Gly - Ser - Gln - Arg - Ala - Leu - Arg - Leu - Ser - Leu - Ala - Gly - Pro - Ile - Leu - Ser - Leu - Cys - Leu - Leu - Pro - Ser - Leu
Species	Human
Reconstitution	For best results and reproducibility, rehydrate the product immediately before use with 50~70 µl DMSO and 25~50 µl distilled water (giving a total volume of 100 µl).
Storage Condition	Upon receipt, store product in its original form at -20°C (or below) and keep away from light, especially after reconstitution.
Content	Each vial contains lyophilized rhodamine labeled rabbit anti-Growth Hormone Secretagogue Receptor type 1b (245-289) (Human) purified IgG.
Recommended Dilution Factor	1:100 (i.e. 10 µl stock in 1 ml dilution buffer) 1:500 (i.e. 2 µl stock in 1 ml dilution buffer)
Absorption	Excitation: 567 nm
Fluorescent_Dye	Rhodamine B

Thursday 20 December, 2012

11985386 requests since Wednesday 05 April, 2006

Copyright © 2012 PHOENIX PHARMACEUTICALS, INC.