

# Growth Hormone Secretagogue Receptor type 1A [Cys0] (330-366) (Human) - Rhodamine Labeled Purified IgG

Catalog#	Standard size	Price	Order
FR-G-001-62	100 $\mu$ l	\$450.00	<input type="button" value="ORDER"/>

## Product description

## References

<b>Catalog #</b>	FR-G-001-62
<b>Standard Size</b>	100 $\mu$ l
<b>Sequence</b>	Cys - Tyr - Arg - Val - Ala - Val - Phe - Arg - Leu - Leu - Gly - Phe - Glu - Pro - Phe - Ser - Gln - Arg - Lys - Leu - Ser - Thr - Leu - Lys - Asp - Glu - Ser - Ser - Arg - Ala - Trp - Thr - Glu - Ser - Ser - Ile - Asn - Thr
<b>Species</b>	Human
<b>Storage Condition</b>	Upon receipt, store product in its original form at -20°C (or below) and keep away from light, especially after reconstitution.
<b>Content</b>	Each vial contains 200 $\mu$ l of Rhodamine-labeled rabbit anti-Growth Hormone Secretagogue Receptor type 1A (330-366) (Human) Purified IgG in liquid form.
<b>Absorption</b>	Excitation: 567 nm
<b>Fluorescent_Dye</b>	Rhodamine B

## Image

