

ANP (Rat, Mouse) - Rhodamine Labeled Purified IgG

Catalog#	Standard size	Price	Order
FR-G-005-24	100 µl	\$450.00	<input type="button" value="ORDER"/>

Product description

More information

References

ANP (Atrial Natriuretic Peptide)

Normal Plasma Peptide/Protein Levels

Catalog #	FR-G-005-24
Standard Size	100 µl
Sequence	Ser - Leu - Arg - Arg - Ser - Ser - Cys - Phe - Gly - Gly - Arg - Ile - Asp - Arg - Ile - Gly - Ala - Gln - Ser - Gly - Leu - Gly - Cys - Asn - Ser - Phe - Arg - Tyr
Species	Rat, Mouse
Storage Condition	Upon receipt, store product in its original form at -20°C (or below) and keep away from light, especially after reconstitution.
Content	Each vial contains 200 µl of Rhodamine-labeled rabbit anti- Atrial Natriuretic Peptide (ANP) (Rat, Mouse) purified IgG in liquid form.
Absorption	Excitation: 567 nm
Fluorescent_Dye	Rhodamine B
Image	

