

TRPM8 / CMR-1 (656-680) (Rat) - Rhodamine Labeled Purified IgG

Catalog#	Standard size	Price	Order
FR-G-050-50	100 µl	\$500.00	<input type="button" value="ORDER"/>

Product description

More information

References

Icilin (AG-3-5)

Catalog #	FR-G-050-50
Standard Size	100 µl
Sequence	Cys - Leu - Glu - Leu - Ala - Val - Glu - Ala - Thr - Asp - Gln - His - Phe - Ile - Ala - Gln - Pro - Gly - Val - Gln - Asn - Phe - Leu - Ser - Lys
Species	Rat
Reconstitution	For best results and reproducibility, rehydrate the product immediately before use with 50~70 µl DMSO and 25~50 µl distilled water (giving a total volume of 100 µl).
Storage Condition	Upon receipt, store product in its original form at -20°C (or below) and keep away from light, especially after reconstitution.
Content	Each vial contains lyophilized rhodamine labeled rabbit anti-CMR-1 (656-680) (Rat) purified IgG.
Recommended Dilution Factor	1:100 (i.e. 10 µl stock in 1 ml dilution buffer) 1:500 (i.e. 2 µl stock in 1 ml dilution buffer)
Absorption	Excitation: 567 nm
Fluorescent_Dye	Rhodamine B

Thursday 20 December, 2012

11985501 requests since Wednesday 05 April, 2006

Copyright © 2012 PHOENIX PHARMACEUTICALS, INC.