

LL37 (Human) - Rhodamine Labeled Purified IgG

Catalog#	Standard size	Price	Order
FR-G-075-06	100 µl	\$500.00	<input type="button" value="ORDER"/>

- Product description**
- More information
- References

LL37 , KR-12 and WLBU2

Anti-Microbial Peptides

KR-12, An Antibacterial Peptide

Normal Plasma Peptide/Protein Levels

Catalog #	FR-G-075-06
Standard Size	100 µl
Sequence	Leu - Leu - Gly - Asp - Phe - Phe - Arg - Lys - Ser - Lys - Glu - Lys - Ile - Gly - Lys - Glu - Phe - Lys - Arg - Ile - Val - Gln - Arg - Ile - Lys - Asp - Phe - Leu - Arg - Asn - Leu - Val - Pro - Arg - Thr - Glu - Ser
Species	Human
Storage Condition	Upon receipt, store product in its original form at -20°C (or below) and keep away from light, especially after reconstitution.
Content	Each vial contains 100 µl of Rhodamine-labeled rabbit anti- LL37 (Human) purified IgG in liquid form.
Absorption	Excitation: 567 nm
Fluorescent_Dye	Rhodamine B
Image	

