

Sef (Human) - Rhodamine Labeled Purified IgG

Catalog#	Standard size	Price	Order
FR-G-SEF-01	100 µl	\$500.00	<input type="button" value="ORDER"/>

Product description

References

Catalog #	FR-G-SEF-01
Standard Size	100 µl
Species	Human
Reconstitution	For best results and reproducibility, rehydrate the product immediately before use with 50~70 µl DMSO and 25~50 µl distilled water (giving a total volume of 100 µl).
Storage Condition	Upon receipt, store product in its original form at -20°C (or below) and keep away from light, especially after reconstitution.
Content	Each vial contains lyophilized rhodamine labeled rabbit anti-sef (human) purified IgG.
Absorption	Excitation: 567 nm
Fluorescent_Dye	Rhodamine B

Thursday 20 December, 2012

11985528 requests since Wednesday 05 April, 2006

Copyright © 2012 PHOENIX PHARMACEUTICALS, INC.