

# Augurin, Prepro (133-148) (Human) - Purified IgG Antibody

Catalog#	Standard size	Price	Order
G-012-24	100 µg	\$250.00	<input type="button" value="ORDER"/>

## Product description

## More information

[Augurin - Research Abstracts, Antibody Stainings, Kits Standard Curves, Sequences, etc.](#)

[Isolation and Identification of Augurin Peptide from the Bovine Hypothalamus](#)

**Catalog #** G-012-24


**Standard Size** 100 µg

**Species** Human

**Host** Rabbit

**Storage Condition** For optimal results, use the antibody as soon as possible after reconstitution. Store in lyophilized form unless needed and reconstitute immediately before use. Once reconstituted, the antibody should be stable for a few days at -4°C. For storage up to a few months, prepare small aliquots after reconstitution and freeze at -20°C or -80°C. Repeated freeze thaw cycles should be strictly avoided.

## Image



Lane #1: Porcine Pancreas Extracts

Lane #2: Bovine Hypothalamus Extracts

Lane #3: Bovine Pituitary Extracts

Lane #4: Porcine Liver Extracts

The major band in the porcine pancreas is 14.2 kDa augurin, but the bovine hypothalamus and pituitary show augurin derived peptides with the size at 14.1, 9.6 kDa and 6.5 kDa which can be detected by anti-augurin (133-148) antibody.