

GLP-1 (7-37) (Human, Rat, Mouse, Porcine, Bovine, Canine, Ovine) - Purified IgG Antibody

Catalog#	Standard size	Price	Order
G-028-13	400 µg	\$325.00	<input type="button" value="ORDER"/>

Product description

More information

References

GLP-1 & GLP-2

Normal Plasma Peptide/Protein Levels

Neuron Function Modulation

Catalog #	G-028-13
Standard Size	400 µg
Sequence	His - Ala - Glu - Gly - Thr - Phe - Thr - Ser - Asp - Val - Ser - Ser - Tyr - Leu - Glu - Gly - Gln - Ala - Ala - Lys - Glu - Phe - Ile - Ala - Trp - Leu - Val - Lys - Gly - Arg - Gly
Species	Human, Rat, Mouse, Porcine, Bovine, Ovine, Canine
Host	Rabbit
Storage Condition	For optimal results, use the antibody as soon as possible after reconstitution. Store in lyophilized form unless needed and reconstitute immediately before use. Once reconstituted, the antibody should be stable for a few days at -4°C. For storage up to a few months, prepare small aliquots after reconstitution and freeze at -20°C or -80°C. Repeated freeze thaw cycles should be strictly avoided.
Contents	Each vial contains 400µg of rabbit IgG in 400µl PBS (137mM NaCl, 2.7mM KCl, 4.3mM Na ₂ HPO ₄ , 1.4mM KH ₂ PO ₄ , pH 7.2) in the liquid form.

Cross Reactivity	Peptide	% Cross-reactivity
	GLP-1 (7-37) (Human)	100
	GLP-1 (Human)	100
	GLP-1 (7-36) (Human)	0
	GLP-2 (Human)	0
	Glucagon (Human)	0
	Insulin (Human)	0