

Glicentin-related Polypeptide (21-50) (Human) - Purified IgG Antibody

Catalog#	Standard size	Price	Order
G-028-17	200 µg	\$450.00	<input type="button" value="ORDER"/>

Product description

References

Catalog #	G-028-17
Standard Size	200 µg
Sequence	Arg - Ser - Leu - Gln - Asp - Thr - Glu - Glu - Lys - Ser - Arg - Ser - Phe - Ser - Ala - Ser - Gln - Ala - Asp - Pro - Leu - Ser - Asp - Pro - Asp - Gln - Met - Asn - Glu - Asp
Species	Human
Host	Rabbit
Storage Condition	For optimal results, use the antibody as soon as possible after reconstitution. Store in lyophilized form unless needed and reconstitute immediately before use. Once reconstituted, the antibody should be stable for a few days at -4°C. For storage up to a few months, prepare small aliquots after reconstitution and freeze at -20°C or -80°C. Repeated freeze thaw cycles should be strictly avoided.
Contents	This vial contains 200µg of rabbit IgG in 200µl PBS (137 mM NaCl, 2.7 mM KCl, 4.3 mM Na ₂ HPO ₄ , 1.4 mM KH ₂ PO ₄ , pH7.5) with 0.1% BSA and 0.05% Thimersol in the liquid form.

Thursday 20 December, 2012

11985929 requests since Wednesday 05 April, 2006

Copyright © 2012 PHOENIX PHARMACEUTICALS, INC.