

Exendin-4 (Heloderma suspectum) - Purified IgG Antibody

Catalog#	Standard size	Price	Order
G-070-94	200 µg	\$450.00	<input type="button" value="ORDER"/>

Product description

More information

References

[Exendin-4](#), [Research Abstracts](#), [Kits Standard Curves](#), [Sequences](#), etc.

[Neuron Function Modulation](#)

Catalog #	G-070-94
Standard Size	200 µg
Sequence	His - Gly - Glu - Gly - Thr - Phe - Thr - Ser - Asp - Leu - Ser - Lys - Gln - Met - Glu - Glu - Glu - Ala - Val - Arg - Leu - Phe - Ile - Glu - Trp - Leu - Lys - Asn - Gly - Gly - Pro - Ser - Ser - Gly - Ala - Pro - Pro - Pro - Ser - NH ₂
Host	Rabbit
Storage Condition	For optimal results, use the antibody as soon as possible after reconstitution. Store in lyophilized form unless needed and reconstitute immediately before use. Once reconstituted, the antibody should be stable for a few days at -4°C. For storage up to a few months, prepare small aliquots after reconstitution and freeze at -20°C or -80°C. Repeated freeze thaw cycles should be strictly avoided.
Contents	This vial contains 200µg of rabbit IgG in 200µl PBS (137mM NaCl, 2.7mM KCl, 4.3mM Na ₂ HPO ₄ , 1.4mM KH ₂ PO ₄ , pH 7.2) in the liquid form.

Thursday 20 December, 2012

11986193 requests since Wednesday 05 April, 2006

Copyright © 2012 PHOENIX PHARMACEUTICALS, INC.