

Augurin, Prepro (133-148) (Human) - Antibody for Immunohistochemistry

Catalog#	Standard size	Price	Order
H-012-24	100 μ l	\$250.00	<input type="button" value="ORDER"/>


Product description

More information


[Augurin - Research Abstracts, Antibody Stainings, Kits Standard Curves, Sequences, etc.](#)

[Isolation and Identification of Augurin Peptide from the Bovine Hypothalamus](#)

Image



Human pancreas tissue was stained by rabbit anti-Augurin (133-148) (H) Antiserum (Cat. No.: H-012-24)



Human pancreas tissue was stained by rabbit anti-Augurin, Prepro (133-148) (H) Antiserum (Cat. No.: H-012-24)

Protocol for Immunohistochemistry

Tissue Sample	Human pancreas tissues
Fixative	10% Formalin
Embedding	Paraffin
Negative control	No primary antibody
Pretreatment	N/A
Blocking	3% H ₂ O ₂ , 2% Normal Goat Serum
Primary Antibody	Rabbit Anti-Augurin, Prepro (133-148) (H) Antiserum (Catalog No.: H-012-24)
Optimal Dilution	1:500, 1 hour at RT
Secondary Antibody	Goat anti-Rabbit IgG, Biotinylated (1:400, 30 min)
Amplification	Streptavidin-HRP (Vector), 1:400, 30 min
Detection system	HRP
Substrate	DAB (Sigma), 3 min
Counterstained	Hematoxylin, 30 Sec

freeze at -20°C or -80°C. Repeated freeze thaw cycles should be strictly avoided.

Recommended Dilution Factor Immunohistochemistry - 1:500

Thursday 20 December, 2012

11986433 requests since Wednesday 05 April, 2006

Copyright © 2012 PHOENIX PHARMACEUTICALS, INC.