

# Neuromedin U, prepro (104-136) (Human) - Antibody for Immunohistochemistry

Catalog#	Standard size	Price	Order
H-045-96	100 µl	\$450.00	<input type="button" value="ORDER"/>

## Product description

## More information

### Neuromedin U (NmU)

<b>Catalog #</b>	H-045-96
<b>Standard Size</b>	100 µl
<b>Sequence</b>	Phe - Leu - Phe - His - Tyr - Ser - Lys - Thr - Gln - Lys - Leu - Gly - Lys - Ser - Asn - Val - Val - Ser - Ser - Val - Val - His - Pro - Leu - Leu - Gln - Leu - Val - Pro - His - Leu - His - Glu
<b>Species</b>	Human
<b>Host</b>	Rabbit
<b>Reconstitution</b>	For best and reproducible results, reconstitute with 100µl of distilled water for the equivalent of undiluted antiserum, immediately before use. Do not refreeze any unused portions.
<b>Storage Condition</b>	For optimal results, use the antibody as soon as possible after reconstitution. Store in lyophilized form unless needed and reconstitute immediately before use. Once reconstituted, the antibody should be stable for a few days at -4°C. For storage up to a few months, prepare small aliquots after reconstitution and freeze at -20°C or -80°C. Repeated freeze thaw cycles should be strictly avoided.
<b>Content</b>	This vial contains 100µl of Rabbit Anti -Prepro-Neuromedin U (104-136) / NmU (104-136) (Human) Serum in the lyophilized form. Prepro-Neuromedin U (104-136) / NmU (104-136) (Human) can be detected with this IgG (1:6400, 0.156µg/ml) in ELISA.

Thursday 20 December, 2012

11986623 requests since Wednesday 05 April, 2006

Copyright © 2012 PHOENIX PHARMACEUTICALS, INC.