

# Neuromedin U-25 (Human) - Antibody for Immunohistochemistry

Catalog#	Standard size	Price	Order
H-046-42	50 µl	\$450.00	<input type="button" value="ORDER"/>

## Product description

## More information

### Neuromedin U (NmU)

<b>Catalog #</b>	H-046-42
<b>Standard Size</b>	50 µl
<b>Sequence</b>	Phe - Arg - Val - Asp - Glu - Glu - Phe - Gln - Ser - Pro - Phe - Ala - Ser - Gln - Ser - Arg - Gly - Tyr - Phe - Leu - Phe - Arg - Pro - Arg - Asn - NH <sub>2</sub>
<b>Species</b>	Human
<b>Host</b>	Rabbit
<b>Reconstitution</b>	For best and reproducible results, reconstitute with 50µl of distilled water for the equivalent of undiluted antiserum, immediately before use. Do not refreeze any unused portions.
<b>Storage Condition</b>	For optimal results, use the antibody as soon as possible after reconstitution. Store in lyophilized form unless needed and reconstitute immediately before use. Once reconstituted, the antibody should be stable for a few days at -4°C. For storage up to a few months, prepare small aliquots after reconstitution and freeze at -20°C or -80°C. Repeated freeze thaw cycles should be strictly avoided.
<b>Content</b>	This vial contains 50µl of Rabbit Anti- Neuromedin U-25 (Human) Serum in the lyophilized form.
<b>Recommended Dilution Factor</b>	Immunofluorescence: 1:200 PAP or ABC: 1:500

Thursday 20 December, 2012

11986629 requests since Wednesday 05 April, 2006

Copyright © 2012 PHOENIX PHARMACEUTICALS, INC.