

Exendin-4 (Heloderma suspectum) - Antibody for Immunohistochemistry

Catalog#	Standard size	Price	Order
H-070-94	100 µl	\$450.00	<input type="button" value="ORDER"/>

Product description

More information

[Exendin-4](#), [Research Abstracts](#), [Kits](#) [Standard Curves](#), [Sequences](#), etc.

IL-6

[Neuron Function Modulation](#)

Catalog #	H-070-94
Standard Size	100 µl
Sequence	His - Gly - Glu - Gly - Thr - Phe - Thr - Ser - Asp - Leu - Ser - Lys - Gln - Met - Glu - Glu - Glu - Ala - Val - Arg - Leu - Phe - Ile - Glu - Trp - Leu - Lys - Asn - Gly - Gly - Pro - Ser - Ser - Gly - Ala - Pro - Pro - Pro - Ser - NH ₂
Specific Specie	Heloderma suspectum
Storage Condition	For optimal results, use the antibody as soon as possible after reconstitution. Store in lyophilized form unless needed and reconstitute immediately before use. Once reconstituted, the antibody should be stable for a few days at -4°C. For storage up to a few months, prepare small aliquots after reconstitution and freeze at -20°C or -80°C. Repeated freeze thaw cycles should be strictly avoided.

Thursday 20 December, 2012

11986724 requests since Wednesday 05 April, 2006

Copyright © 2012 PHOENIX PHARMACEUTICALS, INC.