

Lavender Vacutainer tubes - 100 per Box

Catalog#	Standard size	Price	Order
VT-6450	100 Tubes per Box	\$72.00	<input type="button" value="ORDER"/>

Product description

More information

[Suggested Methods of Peptide Extraction from Tissues in Immunoassay Studies](#)

[Suggested General Method for Blood Collection](#)

[Plasma Extraction](#)

[CSF Extraction](#)

Sample Preparation	Blood collection Plasma extraction Tissue extraction CSF extraction
Catalog #	VT-6450
Standard Size	13 x 100 mm
Radioactivity	This kit contains 1.5 μ Ci of Iodine 125
Range	Maximum Draw Volume: 6 ml

Thursday 20 December, 2012

11986814 requests since Wednesday 05 April, 2006

Copyright © 2012 PHOENIX PHARMACEUTICALS, INC.