

***Rat Anti-Mouse B220 (B cell specific isoform of CD45)  
(rAmB220)***

**Certificate of Analysis and Data Sheet**

<b><i>CLONE</i></b>	mB220
<b><i>IMMUNOGEN</i></b>	Purified mouse LN B cells
<b>Ig SUBCLASS</b>	Rat IgM
<b>PURIFICATION METHOD</b>	Ion exchange column
<b>PROTEIN CONCENTRATION</b>	1 mg/ml in PBS (after reconstitution)
<b>STORAGE PROCEDURES</b>	Lyophilized: store at 4°C. After reconstitution, if not intended for use within a month, aliquot and store at 20°C.
<b>RECONSTITUTION</b>	Reconstitute with of H <sub>2</sub> O. Mix gently, wash the sides of the vial and wait 30-60 seconds before use.
<b>STABILITY / SHELF LIFE</b>	Two years lyophilized, one month in solution at 4°C.
<b>ANTIBODY PROPERTIES</b>	Staining antibody. For staining, use 10 ml/10 <sup>6</sup> cells.
<b>AVAILABLE CONJUGATES</b>	This antibody is only available as purified antibody.