

Phospholipid-TGT

3 mL

0.5 mM Phospholipid Emulsion

For In Vitro Research Use Only

REF PL604T



2022 10

LOT PLS0916



2-8 °C, Do not freeze



Swirl gently before use

Composition

Phospholipid emulsion containing a mixture of highly purified phosphatidyl choline (PC), phosphatidyl serine (PS) and sphingomyelin (SM).

Molar % - PC 42%: PS 28%: SM 30%

Medium: 0.05 mol/L Tris-HCl, pH 7.6.

Properties

A long term stabilized phospholipid emulsion with high procoagulant activity.

Stability

A sealed vial is stable at 2-8°C until the expiry date printed on the label. An opened vial is stable for one month at 2-8°C, provided that any microbial contamination is avoided during use. Do not freeze.

Application

This phospholipid emulsion is suitable for use in hemostasis research related to procoagulant and anticoagulant pathways such as in thrombin generation assays with and without activated Protein C.

Package

10 x 3 mL, Swirl gently before use.



Rossix AB
Taljegardsgatan 3B
SE-431 53 Molndal
Sweden

Phone: +46 31 706 8965 8965
Fax: +46 31 706 8966
info@rossix.com
www.rossix.com