

**Product Name:** 8-Hydroxy-PIPAT oxalate

**Catalog No.:** 0797

**Batch No.:** 2

CAS Number: 1451210-48-2

IUPAC Name: (RS)-trans-8-Hydroxy-2-[N-n-propyl-N-(3'-iodo-2'-propenyl)amino]tetralin oxalate

**1. PHYSICAL AND CHEMICAL PROPERTIES**

**Batch Molecular Formula:** C<sub>16</sub>H<sub>22</sub>INO.C<sub>2</sub>H<sub>2</sub>O<sub>4</sub>


**Batch Molecular Weight:** 461.3

**Physical Appearance:** Off-white solid

**Solubility:** ethanol to 100 mM  
DMSO to 100 mM

**Storage:** Desiccate at -20°C

**Batch Molecular Structure:**



**2. ANALYTICAL DATA**

**TLC:** R<sub>f</sub> = 0.45 (Ethyl acetate:Petroleum ether [1:2])

**Melting Point:** Greater than 75°C

**HPLC:** Shows >98.1% purity

**<sup>1</sup>H NMR:** Consistent with structure

**Microanalysis:**

	Carbon	Hydrogen	Nitrogen
Theoretical	46.87	5.24	3.04
Found	46.76	5.42	3.01

Caution - Not Fully Tested • Research Use Only • Not For Human or Veterinary Use

**Product Name:** 8-Hydroxy-PIPAT oxalate

**Catalog No.:** 0797

**Batch No.:** 2

CAS Number: 1451210-48-2

IUPAC Name: (RS)-trans-8-Hydroxy-2-[N-n-propyl-N-(3'-iodo-2'-propenyl)amino]tetralin oxalate

**Description:**

High affinity 5-HT<sub>1A</sub> receptor agonist ( $K_d = 0.38$  nM). 7-Hydroxy-PIPAT maleate (Cat. No. 0719) also available.

**Physical and Chemical Properties:**


Batch Molecular Formula: C<sub>16</sub>H<sub>22</sub>INO.C<sub>2</sub>H<sub>2</sub>O<sub>4</sub>

Batch Molecular Weight: 461.3

Physical Appearance: Off-white solid

**Minimum Purity:** >98%

**Batch Molecular Structure:**



**Storage:** Desiccate at -20°C

CAUTION - This product is light sensitive and we recommend that the solid material and any solutions obtained are protected from exposure to light.

**Solubility & Usage Info:**

ethanol to 100 mM

DMSO to 100 mM

**Stability and Solubility Advice:**

Some solutions can be difficult to obtain and can be encouraged by rapid stirring, sonication or gentle warming (in a 45-60°C water bath).

Information concerning product stability, particularly in solution, has rarely been reported and in most cases we can only offer a general guide. Our standard recommendations are:

**SOLIDS:** Provided storage is as stated on the product label and the vial is kept tightly sealed, the product can be stored for up to 6 months from date of receipt.

**SOLUTIONS:** We recommend that stock solutions, once prepared, are stored aliquoted in tightly sealed vials at -20°C or below and used within 1 month. Wherever possible solutions should be made up and used on the same day.

**References:**

Zhuang *et al* (1993) Synthesis of (R,S)-trans-8-Hydroxy-2-[N-n-propyl-N-(3'-iodo-2'-propenyl)amino]tetralin (trans-8-OH-PIPAT): a new 5-HT<sub>1A</sub> receptor ligand. J.Med.Chem. **36** 3161. PMID: 8230102.

Caution - Not Fully Tested • Research Use Only • Not For Human or Veterinary Use

Tocris Bioscience is an R&D Systems company

USA & CANADA Tel: (800) 343-7475 EUROPE Tel: +44 (0)1235 529449 CHINA Tel: +86 (21) 52380373

[www.RnDSystems.com](http://www.RnDSystems.com)