

Product Name: Ketanserin tartrate

Catalog No.: 0908

Batch No.: 4

CAS Number: 83846-83-7

EC Number: 281-062-8

IUPAC Name: 3-[2-[4-(4-Fluorobenzoyl)-1-piperidinyl]ethyl]-2,4[1*H*,3*H*]-quinazolin-2(1*H*)-one tartrate

1. PHYSICAL AND CHEMICAL PROPERTIES

Batch Molecular Formula: $C_{22}H_{22}FN_3O_3 \cdot C_4H_6O_6 \cdot H_2O$


Batch Molecular Weight: 563.54

Physical Appearance: White solid

Solubility: water to 10 mM
DMSO to 100 mM

Storage: Store at RT

Batch Molecular Structure:



2. ANALYTICAL DATA

Melting Point: Between 187 - 189°C

HPLC: Shows 99% purity

¹H NMR: Consistent with structure

Mass Spectrum: Consistent with structure

Microanalysis:

	Carbon	Hydrogen	Nitrogen
Theoretical	55.42	5.37	7.46
Found	55.35	5.53	7.46

Caution - Not Fully Tested • Research Use Only • Not For Human or Veterinary Use

Product Name: Ketanserin tartrate

Catalog No.: 0908

Batch No.: 4

CAS Number: 83846-83-7

EC Number: 281-062-8

IUPAC Name: 3-[2-[4-(4-Fluorobenzoyl)-1-piperidinyl]ethyl]-2,4[1*H*,3*H*]-quinazolin-1(2*H*)-one tartrate

Description:

Selective 5-HT_{2A} serotonin receptor antagonist; can also be used to discriminate between 5-HT_{1D} and 5-HT_{1B} receptor subtypes.

Physical and Chemical Properties:


Batch Molecular Formula: C₂₂H₂₂FN₃O₃·C₄H₆O₆·H₂O

Batch Molecular Weight: 563.54

Physical Appearance: White solid

Minimum Purity: >97%

Batch Molecular Structure:



Storage: Store at RT

Solubility & Usage Info:

water to 10 mM
DMSO to 100 mM

Stability and Solubility Advice:

Some solutions can be difficult to obtain and can be encouraged by rapid stirring, sonication or gentle warming (in a 45-60°C water bath).

Information concerning product stability, particularly in solution, has rarely been reported and in most cases we can only offer a general guide. Our standard recommendations are:

SOLIDS: Provided storage is as stated on the product label and the vial is kept tightly sealed, the product can be stored for up to 6 months from date of receipt.

SOLUTIONS: We recommend that stock solutions, once prepared, are stored aliquoted in tightly sealed vials at -20°C or below and used within 1 month. Wherever possible solutions should be made up and used on the same day.

References:

Leyson *et al* (1981) Receptor binding profile of R41 468, a novel antagonist at 5-HT₂ receptors. *Life Sci.* **28** 1015. PMID: 6261070.

Van Nueten *et al* (1981) Vascular effects of ketanserin (R41 468), a novel antagonist of 5-HT₂ serotonergic receptors. *J.Pharmacol.Exp.Ther.* **218** 217. PMID: 6113280.

Zgombick *et al* (1995) Ketanserin and ritanserin discriminate between recombinant human 5-HT_{1Dα} and 5-HT_{1Dβ} receptor subtypes. *Eur.J.Pharmacol.* **291** 9. PMID: 8549648.

Razzaque *et al* (1995) Differences in the effects of ketanserin and GR 127935 on 5-HT-receptor mediated responses in rabbit saphenous vein and guinea-pig jugular vein. *Eur.J.Pharmacol.* **283** 199. PMID: 7498311.

Caution - Not Fully Tested • Research Use Only • Not For Human or Veterinary Use

Tocris Bioscience is an R&D Systems company
USA & CANADA Tel: (800) 343-7475 EUROPE Tel: +44 (0)1235 529449 CHINA Tel: +86 (21) 52380373
www.RnDSystems.com