

**Product Name:** CGRP 8-37 (rat)

**Catalog No.:** 1169

**Batch No.:** 24

**CAS Number:** 129121-73-9

**1. PHYSICAL AND CHEMICAL PROPERTIES**

**Batch Molecular Formula:** C<sub>138</sub>H<sub>224</sub>N<sub>42</sub>O<sub>41</sub>  
**Batch Molecular Weight:** 3127.5  
**Physical Appearance:** White lyophilised solid  
**Net Peptide Content:** 82.9%  
**Counter Ion:** TFA  
**Solubility:** Soluble to 1 mg/ml in water  
**Storage:** Desiccate at -20°C  
**Peptide Sequence:** Val-Thr-His-Arg-Leu-Ala-Gly-Leu-Leu-Ser-Arg-Ser-Gly-Gly-Val-Val-Lys-Asp-Asn-Phe-Val-Pro-Thr-Asn-Val-Gly-Ser-Glu-Ala-Phe-NH<sub>2</sub>

**2. ANALYTICAL DATA**

**HPLC:** Shows 96.3% purity  
**Mass Spectrum:** Consistent with structure

**3. AMINO ACID ANALYSIS DATA**

Amino Acid		Theoretical	Actual	Amino Acid		Theoretical	Actual
Ala	2.00	1.90	Lys	1.00	1.07		
Arg	2.00	2.08	Met				
Asx	3.00	3.11	Phe	2.00	2.02		
Cys			Pro	1.00	1.03		
Glx	1.00	1.01	Ser	3.00	2.70		
Gly	4.00	4.00	Thr	2.00	1.81		
His	1.00	1.03	Trp				
Ile			Tyr				
Leu	3.00	3.13	Val	5.00	4.86		

Caution - Not Fully Tested • Research Use Only • Not For Human or Veterinary Use

**Product Name:** CGRP 8-37 (rat)

**Catalog No.:** 1169

**Batch No.:** 24

CAS Number: 129121-73-9

**Description:**

Peptide antagonist for CGRP<sub>1</sub> receptors.

**Physical and Chemical Properties:**

Batch Molecular Formula: C<sub>138</sub>H<sub>224</sub>N<sub>42</sub>O<sub>41</sub>

Batch Molecular Weight: 3127.5

Physical Appearance: White lyophilised solid

**Peptide Sequence:**

Val-Thr-His-Arg-Leu-Ala-Gly-Leu-Leu-Ser-  
Arg-Ser-Gly-Gly-Val-Val-Lys-Asp-Asn-Phe-  
Val-Pro-Thr-Asn-Val-Gly-Ser-Glu-Ala-Phe-NH<sub>2</sub>

**Storage:** Desiccate at -20°C

**Solubility & Usage Info:**

Soluble to 1 mg/ml in water

This product is supplied as a lyophilized solid and may be very hard to visualize. Solutions should be made by adding solvent directly to the vial. The vial should then be vortexed vigorously to ensure the product has completely dissolved.

**Net Peptide Content:** 82.9% (Remaining weight made up of counterions and residual water).

**Counter Ion:** TFA

**Stability and Solubility Advice:**

Some solutions can be difficult to obtain and can be encouraged by rapid stirring, sonication or gentle warming (in a 45-60°C water bath).

Peptides in solution are much less stable than in lyophilized form. This is especially true for peptides whose sequences contain amino acids such as Cys, Met, Trp, Asn, Gln, and N-terminal Glu.

Therefore we recommend storing peptides in solution for as short a time as possible. Avoid repeated freeze thaw cycles by dividing the peptide solution into aliquots and storing the aliquots at -20°C. Any portion of an aliquot unused after thawing should be discarded.

Peptides stored in solution can occasionally be susceptible to bacterial degradation. We recommend using sterile solutions or passing the peptide solution through a 0.2 µm filter to remove potential bacterial contamination whenever possible.

**References:**

**Poyner** (1995) Pharmacology of receptors for calcitonin gene-related peptide and amylin. *TIPS* **16** 424. PMID: 8578616.

**Poyner et al** (1998) Structural determinants for binding to CGRP receptors expressed by human SK-N-MC and Col 29 cells: studies with chimeric and other peptides. *Br.J.Pharmacol.* **124** 1659. PMID: 9756381.

**Wisskirchen et al** (1998) Pharmacological characterization of CGRP receptors mediating relaxation of the rat pulmonary artery and inhibition of twitch responses of the rat vas deferens. *Br.J.Pharmacol.* **123** 1673. PMID: 9605575.

Caution - Not Fully Tested • Research Use Only • Not For Human or Veterinary Use

Tocris Bioscience is an R&D Systems company

USA & CANADA Tel: (800) 343-7475 EUROPE Tel: +44 (0)1235 529449 CHINA Tel: +86 (21) 52380373

www.RnDSystems.com

