

Product Name: PHM 27 (human)

Catalog No.: 2080

Batch No.: 1

CAS Number: 118025-43-7

1. PHYSICAL AND CHEMICAL PROPERTIES

Batch Molecular Formula: C₁₃₅H₂₁₄N₃₄O₄₀S
Batch Molecular Weight: 2985.44
Physical Appearance: White lyophilised solid
Net Peptide Content: 88%
Solubility: Soluble to 1 mg/ml in 5% acetonitrile / water
Storage: Desiccate at -20°C
Peptide Sequence: His-Ala-Asp-Gly-Val-Phe-Thr-Ser-Asp-Phe-
Ser-Lys-Leu-Leu-Gly-Gln-Leu-Ser-Ala-Lys-
Lys-Tyr-Leu-Glu-Ser-Leu-Met-NH₂

2. ANALYTICAL DATA

HPLC: Shows >95% purity
Mass Spectrum: Consistent with structure

Caution - Not Fully Tested • Research Use Only • Not For Human or Veterinary Use

Tocris Bioscience is an R&D Systems company
USA & CANADA Tel: (800) 343-7475 EUROPE Tel: +44 (0)1235 529449 CHINA Tel: +86 (21) 52380373
www.RnDSystems.com

R&D
SYSTEMS®

Product Name: PHM 27 (human)

Catalog No.: 2080

Batch No.: 1

CAS Number: 118025-43-7

Description:

Endogenous peptide product of human prepro-VIP and analog of porcine PHI-27; potent agonist for the human calcitonin receptor (EC_{50} = 11 nM). Transgenic mice expressing the human VIP/PHM 27 gene in pancreatic β -islets display enhanced glucose-induced insulin secretion.

Physical and Chemical Properties:

Batch Molecular Formula: $C_{135}H_{214}N_{34}O_{40}S$

Batch Molecular Weight: 2985.44

Physical Appearance: White lyophilised solid

Peptide Sequence:

His-Ala-Asp-Gly-Val-Phe-Thr-Ser-Asp-Phe-
Ser-Lys-Leu-Leu-Gly-Gln-Leu-Ser-Ala-Lys-
Lys-Tyr-Leu-Glu-Ser-Leu-Met-NH₂

Storage: Desiccate at -20°C

Solubility & Usage Info:

Soluble to 1 mg/ml in 5% acetonitrile / water

This product is supplied as a lyophilized solid and may be very hard to visualize. Solutions should be made by adding solvent directly to the vial. The vial should then be vortexed vigorously to ensure the product has completely dissolved.

Net Peptide Content: 88% (Remaining weight made up of counterions and residual water).

Stability and Solubility Advice:

Some solutions can be difficult to obtain and can be encouraged by rapid stirring, sonication or gentle warming (in a 45-60°C water bath).

Peptides in solution are much less stable than in lyophilized form. This is especially true for peptides whose sequences contain amino acids such as Cys, Met, Trp, Asn, Gln, and N-terminal Glu.

Therefore we recommend storing peptides in solution for as short a time as possible. Avoid repeated freeze thaw cycles by dividing the peptide solution into aliquots and storing the aliquots at -20°C. Any portion of an aliquot unused after thawing should be discarded.

Peptides stored in solution can occasionally be susceptible to bacterial degradation. We recommend using sterile solutions or passing the peptide solution through a 0.2 μ m filter to remove potential bacterial contamination whenever possible.

References:

Itoh et al (1983) Human preprovasoactive intestinal polypeptide contains a novel PHI-27-like peptide, PHM-27. *Nature* **304** 547. PMID: 6571696.

Kato et al (1994) Transgenic mice overexpressing human vasoactive intestinal peptide (VIP) gene in pancreatic β cells. *J.Biol.Chem.* **269** 21223. PMID: 8063743.

Ma et al (2004) Discovery of novel peptide/receptor interactions: identification of PHM-27 as a potent agonist of the human calcitonin receptor. *Biochem.Pharmacol.* **67** 1279. PMID: 15013843.

Caution - Not Fully Tested • Research Use Only • Not For Human or Veterinary Use

Tocris Bioscience is an R&D Systems company

USA & CANADA Tel: (800) 343-7475 EUROPE Tel: +44 (0)1235 529449 CHINA Tel: +86 (21) 52380373

www.RnDSystems.com

