

Product Name: Leukotriene B₄

Catalog No.: 2307


Batch No.: 8

CAS Number: 71160-24-2

IUPAC Name: (5S,6Z,8E,10E,12R,14Z)-5,12-Dihydroxy-6,8,10,14-eicosatetraenoic acid

1. PHYSICAL AND CHEMICAL PROPERTIES

Batch Molecular Formula: C₂₀H₃₂O₄
Batch Molecular Weight: 336.47
Physical Appearance: Colourless liquid
Solubility: Soluble in ethanol (supplied pre-dissolved in anhydrous ethanol, 50µg/ml)
Storage: Store at -20°C
Batch Molecular Structure:



2. ANALYTICAL DATA

HPLC: Shows 98.8% purity
Mass Spectrum: Consistent with structure

Caution - Not Fully Tested • Research Use Only • Not For Human or Veterinary Use

Product Name: Leukotriene B₄

Catalog No.: 2307

Batch No.: 8

CAS Number: 71160-24-2

IUPAC Name: (5S,6Z,8E,10E,12R,14Z)-5,12-Dihydroxy-6,8,10,14-eicosatetraenoic acid

Description:

Potent lipid inflammatory mediator derived from the 5-lipoxygenase pathway of arachidonic acid metabolism. Binds to BLT₁ and BLT₂ receptors and acts as a potent chemotactic agent and activator of leukocytes. Also displays antiviral activity towards DNA viruses and retroviruses.


Physical and Chemical Properties:

Batch Molecular Formula: C₂₀H₃₂O₄

Batch Molecular Weight: 336.47

Physical Appearance: Colourless liquid

Batch Molecular Structure:



Storage: Store at -20°C

CAUTION - This product is light sensitive and we recommend that the solid material and any solutions obtained are protected from exposure to light.

Solubility & Usage Info:

Soluble in ethanol (supplied pre-dissolved in anhydrous ethanol, 50µg/ml)

CAUTION - This product is supplied dissolved in anhydrous ethanol at a concentration of 50ug/ml. Aqueous solutions should be stored on ice and used within 12 hours of preparation.

Stability and Solubility Advice:

Some solutions can be difficult to obtain and can be encouraged by rapid stirring, sonication or gentle warming (in a 45-60°C water bath).

Information concerning product stability, particularly in solution, has rarely been reported and in most cases we can only offer a general guide. Our standard recommendations are:

SOLIDS: Provided storage is as stated on the product label and the vial is kept tightly sealed, the product can be stored for up to 6 months from date of receipt.

SOLUTIONS: We recommend that stock solutions, once prepared, are stored aliquoted in tightly sealed vials at -20°C or below and used within 1 month. Wherever possible solutions should be made up and used on the same day.

References:

Bray et al (1981) Leukotriene B₄: an inflammatory mediator in vivo. *Prostaglandins* **22** 213. PMID: 6270742.

Ford-Hutchinson (1990) Leukotriene B₄ in inflammation. *Crit.Rev.Immunol.* **10** 1. PMID: 2155000.

Yokomitso et al (2001) Leukotriene B₄: metabolism and signal transduction. *Arch.Biochem.Biophys.* **385** 231. PMID: 11368003.

Caution - Not Fully Tested • Research Use Only • Not For Human or Veterinary Use

Tocris Bioscience is an R&D Systems company

USA & CANADA Tel: (800) 343-7475 EUROPE Tel: +44 (0)1235 529449 CHINA Tel: +86 (21) 52380373

www.RnDSystems.com