

Product Name: Concanamycin A

Catalog No.: 2656


Batch No.: 1

CAS Number: 80890-47-7

IUPAC Name: (3*Z*,5*E*,7*R*,8*R*,9*S*,10*S*,11*R*,13*E*,15*E*,17*S*,18*R*)-18-[(1*S*,2*R*,3*S*)-3-[(2*R*,4*R*,5*S*,6*R*)-4-[[4-*O*-(Aminocarbonyl)-2,6-dideoxy-β-D-arabino-hexopyranosyl]oxy]tetrahydro-2-hydroxy-5-methyl-6-(1*E*)-1-propenyl-2*H*-pyran-2-yl]-2-hydroxy-1-methylbutyl]-9-ethyl-8,10-dihydroxy-3,17-dimethoxy-5,7,11,13-tetramethyloxacyclooctadeca-3,5,13,15-tetraen-2-one

1. PHYSICAL AND CHEMICAL PROPERTIES

Batch Molecular Formula:	C ₄₆ H ₇₅ NO ₁₄
Batch Molecular Weight:	866.09
Physical Appearance:	White lyophilised solid
Solubility:	Soluble in DMSO
Storage:	Store at -20°C
Batch Molecular Structure:	



2. ANALYTICAL DATA

HPLC:	Shows >97% purity
Mass Spectrum:	Consistent with structure

Caution - Not Fully Tested • Research Use Only • Not For Human or Veterinary Use

Product Name: Concanamycin A

Catalog No.: 2656

Batch No.: 1

CAS Number: 80890-47-7

IUPAC Name: (3Z,5E,7R,8R,9S,10S,11R,13E,15E,17S,18R)-18-[(1S,2R,3S)-3-[(2R,4R,5S,6R)-4-[[4-O-(Aminocarbonyl)-2,6-dideoxy-β-D-arabino-hexopyranosyl]oxy]tetrahydro-2-hydroxy-5-methyl-6-(1E)-1-propenyl-2H-pyran-2-yl]-2-hydroxy-1-methylbutyl]-9-ethyl-8,10-dihydroxy-3,17-dimethoxy-5,7,11,13-tetramethyloxacyclooctadeca-3,5,13,15-tetraen-2-one

Description:

Specific inhibitor of V-type (vacuolar) H⁺-ATPase that displays > 2000-fold selectivity over other H⁺-ATPases (IC₅₀ values are 9.2, > 20000, > 20000 nM for yeast V-type, F-type, P-type H⁺-ATPases and porcine P-type Na⁺,K⁺-ATPase respectively). Blocks cell surface expression of virus envelope glycoproteins without affecting synthesis and exhibits cytotoxicity in several cell lines.

Physical and Chemical Properties:


Batch Molecular Formula: C₄₆H₇₅NO₁₄

Batch Molecular Weight: 866.09

Physical Appearance: White lyophilised solid

Minimum Purity: >90%

Batch Molecular Structure:



References:

Muroi et al (1993) Folimycin (Concanamycin A), a specific inhibitor of V-ATPase, blocks intracellular translocation of the glycoprotein of vesicular stomatitis virus before arrival to the golgi apparatus. *Cell Struct.Function.* **18** 139.

Nishihara et al (1995) Specific inhibitors of vacuolar type H⁺-ATPases induce apoptotic cell death. *Biochem.Biophys.Res.Comm.* **212** 255.

Hong et al (2006) Nitric oxide production by the vacuolar-type (H⁺)-ATPase inhibitors bafilomycin A1 and concanamycin A and its possible role in apoptosis in RAW 264.7 cells. *J.Pharmacol.Exp.Ther.* **319** 672. PMID: 16895977.

Storage: Store at -20°C

CAUTION - This product is light sensitive and we recommend that the solid material and any solutions obtained are protected from exposure to light.

Solubility & Usage Info:

Soluble in DMSO

Lyophilized solids can be hard to visualize therefore solutions should be made by adding solvent directly to the vial. The vial should then be vortexed vigorously to ensure the product has completely dissolved.

Stability and Solubility Advice:

Some solutions can be difficult to obtain and can be encouraged by rapid stirring, sonication or gentle warming (in a 45-60°C water bath).

Information concerning product stability, particularly in solution, has rarely been reported and in most cases we can only offer a general guide. Our standard recommendations are:

SOLIDS: Provided storage is as stated on the product label and the vial is kept tightly sealed, the product can be stored for up to 6 months from date of receipt.

SOLUTIONS: We recommend that stock solutions, once prepared, are stored aliquoted in tightly sealed vials at -20°C or below and used within 1 month. Wherever possible solutions should be made up and used on the same day.

Caution - Not Fully Tested • Research Use Only • Not For Human or Veterinary Use

Tocris Bioscience is an R&D Systems company

USA & CANADA Tel: (800) 343-7475 EUROPE Tel: +44 (0)1235 529449 CHINA Tel: +86 (21) 52380373

www.RnDSystems.com