

MOUSE MONOCLONAL ANTIBODY AGAINST PIG INTERFERON ALPHA

Product Number: **27100-1**

Size: **500 mg**



Product Information

Description: Mouse Anti-Pig Interferon Alpha, neutralizing

Clone: K9

Lot Number: #

Quantity Shipped: 1 vial of 500 µg

Concentration: 1.5 mg/ml

Buffer: 0.1 M Acetate, 50% glycerol, pH 7.5

Antigen: pig interferon alpha

Isotype: IgG₁

Bioactivity: neutralizes human pig alpha

Usage: Western blot; immunoprecipitation; immunohistochemistry; ELISA

References: L'Haridon *et al.* (1991) *Hybridoma*, **10**: 35-47.
Splichal *et al.* (1994) *Immunology Letters*, **43**: 203-208. (ELISA)

Shipping & Storage Information

Date Shipped:

Shipping Temperature: Wet ice

Physical State of Product During Shipping: Frozen

Special Conditions/Comments: After receipt, this product should be kept at -20EC for retention of full activity. When thawing, the contents of the tube should be apportioned in separate tubes so that freezing and thawing is kept to a minimum.

This Product Is For Research Use Only. It Is Not For Human Use.

PBL Biomedical Laboratories

Toll Free: 877-PBL-8881

Tel: 732-828-8881

Fax: 732-828-3736

Email: info@interferonsource.com