

# MOUSE MONOCLONAL ANTIBODY AGAINST PIG INTERFERON ALPHA

Product Number: **27105-1**

Size: **500 mg**



## Product Information

*Description:* Mouse Anti-Pig Interferon Alpha, neutralizing

*Clone:* F17

*Lot Number:* #

*Quantity Shipped:* 1 vial of 500 µg

*Concentration:* 1.5 mg/ml

*Buffer:* 0.1 M Acetate, 50% glycerol, pH 7.5

*Antigen:* pig interferon alpha

*Isotype:* IgG<sub>1</sub>

*Bioactivity:* neutralizes human pig alpha

*Usage:* Western blot; immunoprecipitation; immunohistochemistry; ELISA

*References:* L'Haridon *et al.* (1991) *Hybridoma*, **10**: 35-47.  
Splichal *et al.* (1994) *Immunology Letters*, **43**: 203-208. (ELISA)

## Shipping & Storage Information

*Date Shipped:*

*Shipping Temperature:* Wet ice

*Physical State of Product During Shipping:* Frozen

*Special Conditions/Comments:* After receipt, this product should be kept at -20EC for retention of full activity. When thawing, the contents of the tube should be apportioned in separate tubes so that freezing and thawing is kept to a minimum.

**This Product Is For Research Use Only. It Is Not For Human Use.**

PBL Biomedical Laboratories

Toll Free: 877-PBL-8881

Tel: 732-828-8881

Fax: 732-828-3736

Email: [info@interferonsource.com](mailto:info@interferonsource.com)