

Product Name: 16,16-Dimethyl Prostaglandin E2

Catalog No.: 4027


Batch No.: 4

CAS Number: 39746-25-3

IUPAC Name: (5Z,11 α ,13E,15R)-11,15-Dihydroxy-16,16-dimethyl-9-oxo-prosta-5,13-dien-1oic acid

1. PHYSICAL AND CHEMICAL PROPERTIES

Batch Molecular Formula: C₂₂H₃₆O₅
Batch Molecular Weight: 380.52
Physical Appearance: Colourless liquid
Solubility: Soluble in methyl acetate (supplied pre-dissolved -10mg/ml)
Storage: Store at -80°C
Batch Molecular Structure:



2. ANALYTICAL DATA

HPLC: Shows 98.1% purity
Mass Spectrum: Consistent with structure

Caution - Not Fully Tested • Research Use Only • Not For Human or Veterinary Use

Product Name: 16,16-Dimethyl Prostaglandin E₂

Catalog No.: 4027

Batch No.: 4

CAS Number: 39746-25-3

IUPAC Name: (5Z,11 α ,13E,15R)-11,15-Dihydroxy-16,16-dimethyl-9-oxo-prosta-5,13-dien-1oic acid

Description:

Synthetic derivative of prostaglandin E₂ (Cat. No. 2296). Increases embryonic stem cell (ESC) hematopoietic colony formation in mouse bone marrow. Also induces an increase in hematopoietic stem cell (HSC) numbers and enhances Wnt activity within the HSC population in zebrafish embryos.

Physical and Chemical Properties:


Batch Molecular Formula: C₂₂H₃₆O₅

Batch Molecular Weight: 380.52

Physical Appearance: Colourless liquid

Minimum Purity: >98%

Batch Molecular Structure:



Storage: Store at -80°C

CAUTION - This product is light sensitive and we recommend that the solid material and any solutions obtained are protected from exposure to light.

Solubility & Usage Info:

Soluble in methyl acetate (supplied pre-dissolved -10mg/ml)

This compound is supplied pre-dissolved in Methyl acetate (10mg/ml). To change the solvent, evaporate the methyl acetate under a gentle stream of nitrogen and immediately add the chosen solvent (preferably purged with nitrogen beforehand). The solubility of 16,16-Dimethyl-PGE₂ is greater than 100mM in both DMSO and Ethanol, and about 10mM in PBS. Solutions in DMSO and Ethanol can be kept at -80°C for at least 6 months. Use aqueous solutions within 12 hours.

Stability and Solubility Advice:

Some solutions can be difficult to obtain and can be encouraged by rapid stirring, sonication or gentle warming (in a 45-60°C water bath).

Information concerning product stability, particularly in solution, has rarely been reported and in most cases we can only offer a general guide. Our standard recommendations are:

SOLIDS: Provided storage is as stated on the product label and the vial is kept tightly sealed, the product can be stored for up to 6 months from date of receipt.

SOLUTIONS: We recommend that stock solutions, once prepared, are stored aliquoted in tightly sealed vials at -20°C or below and used within 1 month. Wherever possible solutions should be made up and used on the same day.

References:

Lord et al (2007) Prostaglandin E₂: making more of your marrow. *Cell Cycle* **6** 3054. PMID: 18075310.

North et al (2007) Prostaglandin E₂ regulates vertebrate haematopoietic stem cell homeostasis. *Cell* **147** 1007. PMID: 17581586.

Goessling et al (2009) Genetic interaction of PGE₂ and wnt signaling regulates developmental specification of stem cells and regeneration. *Cell* **136** 1136. PMID: 19303855.

Caution - Not Fully Tested • Research Use Only • Not For Human or Veterinary Use

Tocris Bioscience is an R&D Systems company

USA & CANADA Tel: (800) 343-7475 EUROPE Tel: +44 (0)1235 529449 CHINA Tel: +86 (21) 52380373

www.RnDSystems.com