


Product Name: 3,5-DHBA
CAS Number: 99-10-5
IUPAC Name: 3,5-Dihydroxybenzoic acid

Catalog No.: 4599 **Batch No.:** 1
EC Number: 202-730-7

1. PHYSICAL AND CHEMICAL PROPERTIES

Batch Molecular Formula: C₇H₆O₄
Batch Molecular Weight: 154.12
Physical Appearance: White solid
Solubility: water to 100 mM
DMSO to 100 mM
ethanol to 100 mM
Storage: Store at RT
Batch Molecular Structure:



2. ANALYTICAL DATA

HPLC: Shows 100% purity
¹H NMR: Consistent with structure
Mass Spectrum: Consistent with structure
Microanalysis:

	Carbon	Hydrogen	Nitrogen
Theoretical	54.55	3.92	
Found	54.4	4.01	

Caution - Not Fully Tested • Research Use Only • Not For Human or Veterinary Use

Product Name: 3,5-DHBA

Catalog No.: 4599

Batch No.: 1

CAS Number: 99-10-5

EC Number: 202-730-7

IUPAC Name: 3,5-Dihydroxybenzoic acid

Description:

Selective agonist of hydroxycarboxylic acid receptor 1 (HCA₁, also known as GPR81) (EC₅₀ ~150 μM). Inhibits lipolysis in wild-type mouse adipocytes. Displays minimal activity at HCA₂ (at concentrations >10 mM) and no activity at HCA₃, free fatty acid receptor (FFAR) 2 or FFAR3.

Physical and Chemical Properties:


Batch Molecular Formula: C₇H₆O₄

Batch Molecular Weight: 154.12

Physical Appearance: White solid

Minimum Purity: >99%

Batch Molecular Structure:



Storage: Store at RT

Solubility & Usage Info:

water to 100 mM
DMSO to 100 mM
ethanol to 100 mM

Stability and Solubility Advice:

Some solutions can be difficult to obtain and can be encouraged by rapid stirring, sonication or gentle warming (in a 45-60°C water bath).

Information concerning product stability, particularly in solution, has rarely been reported and in most cases we can only offer a general guide. Our standard recommendations are:

SOLIDS: Provided storage is as stated on the product label and the vial is kept tightly sealed, the product can be stored for up to 6 months from date of receipt.

SOLUTIONS: We recommend that stock solutions, once prepared, are stored aliquoted in tightly sealed vials at -20°C or below and used within 1 month. Wherever possible solutions should be made up and used on the same day.

References:

Liu *et al* (2012) 3,5-dihydroxybenzoic acid, a specific agonist for hydroxycarboxylic acid 1, inhibits lipolysis in adipocytes. *J.Pharmacol.Exp.Ther.* **341** 794. PMID: 22434674.

Caution - Not Fully Tested • Research Use Only • Not For Human or Veterinary Use

Tocris Bioscience is an R&D Systems company
USA & CANADA Tel: (800) 343-7475 EUROPE Tel: +44 (0)1235 529449 CHINA Tel: +86 (21) 52380373
www.RnDSystems.com

R&D
SYSTEMS®