



Certificate of Analysis

www.tocris.com

Product Name: GsMTx4 Catalog No.: 4912 Batch No.: 4

1. PHYSICAL AND CHEMICAL PROPERTIES

Batch Molecular Formula: $C_{185}H_{273}N_{49}O_{45}S_6$

Batch Molecular Weight: 4095.86

Physical Appearance: White lyophilised solid

Counter Ion: TFA

Solubility: Soluble to 1 mg/ml in water

Storage: Store at -20°C

Peptide Sequence:

Gly-Cys-Leu-Glu-Phe-Trp-Trp-Lys-Cys-Asn-

Pro-Asn-Asp-Asp-Lys-Cys-Cys-Arg-Pro-Lys-

Leu-Lys-Cys-Ser-Lys-Leu-Phe-Lys-Leu-Cys

Asn-Phe-Ser-Phe-NH₂

2. ANALYTICAL DATA

HPLC: Shows >99% purity

Mass Spectrum: Consistent with structure



Product Information

Print Date: Apr 28th 2015

www.tocris.com

Product Name: GsMTx4 Catalog No.: 4912 Batch No.: 4

Description:

Inhibitor of the mechanosensitive TRP channels TRPC1 and TRPC6. Also blocks stretch-activated cation channels in astrocytes, cardiac cells, and smooth and skeletal muscle cells.


Physical and Chemical Properties:

Batch Molecular Formula: C₁₈₅H₂₇₃N₄₉O₄₅S₆

Batch Molecular Weight: 4095.86

Physical Appearance: White lyophilised solid

Peptide Sequence:



Storage: Store at -20°C

Solubility & Usage Info:

Soluble to 1 mg/ml in water

This product is supplied as a lyophilized solid and may be very hard to visualize. Solutions should be made by adding solvent directly to the vial. The vial should then be vortexed vigorously to ensure the product has completely dissolved.

Counter Ion: TFA

Stability and Solubility Advice:

Some solutions can be difficult to obtain and can be encouraged by rapid stirring, sonication or gentle warming (in a 45-60°C water bath).

Peptides in solution are much less stable than in lyophilized form. This is especially true for peptides whose sequences contain amino acids such Cys, Met,Trp, Asn, Gln, and N-terminal Glu.

Therefore we recommend storing peptides in solution for as short a time as possible. Avoid repeated freeze thaw cycles by dividing the peptide solution into aliquots and storing the aliquots at -20°C. Any portion of an aliquot unused after thawing should be discarded.

Peptides stored in solution can occasionally be susceptible to bacterial degradation. We recommend using sterile solutions or passing the peptide solution through a 0.2 μ m filter to remove potential bacterial contamination whenever possible.

References:

Suchyna et al (2000) Identification of a peptide toxin from *Grammostola spatulata* spider venom that blocks cation-selective stretch-activated channels. J.Gen.Physiol. 115 583. PMID: 10779316.

Drew *et al* (2007) High-threshold mechanosensitive ion channels blocked by a novel conopeptide mediate pressure-evoked pain. PLoS One **2** e515. PMID: 17565368.

Bowman *et al* (2007) Mechanosensitive ion channels and the peptide inhibitor GsMTx-4: history, properties, mechanisms and pharmacology. Toxicon. *49* 249. PMID: 17157345.

