

DESCRIPTION

Species Reactivity	Mouse
Specificity	Detects mouse Dkk-1 in ELISAs and Western blots. In sandwich ELISAs, less than 6% cross-reactivity with recombinant human Dkk-1 and less than 0.1% cross-reactivity with recombinant mouse (rm) Dkk-2, rmDkk-3, and rmDkk-4 is observed.
Source	Polyclonal Goat IgG
Purification	Antigen Affinity-purified
Immunogen	Mouse myeloma cell line NS0-derived recombinant mouse Dkk-1 Ser30-His272 Accession # O54908
Formulation	Lyophilized from a 0.2 µm filtered solution in PBS with Trehalose. See Certificate of Analysis for details. *Small pack size (-SP) is supplied as a 0.2 µm filtered solution in PBS.

APPLICATIONS

Please Note: Optimal dilutions should be determined by each laboratory for each application. *General Protocols* are available in the *Technical Information* section on our website.

	Recommended Concentration	Sample
Western Blot	0.1 µg/mL	Recombinant Mouse Dkk-1 (Catalog # 1765-DK)
Immunohistochemistry	5-15 µg/mL	See Below
Mouse Dkk-1 Sandwich Immunoassay		Reagent
ELISA Capture	0.2-0.8 µg/mL	Mouse Dkk-1 Antibody (Catalog # AF1765)
ELISA Detection	0.1-0.4 µg/mL	Mouse Dkk-1 Biotinylated Antibody (Catalog # BAF1765)
Standard		Recombinant Mouse Dkk-1 (Histidine-tagged) (Catalog # 1765-DK)

DATA

Immunohistochemistry



Dkk-1 in Mouse Embryo. Dkk-1 was detected in immersion fixed frozen sections of mouse embryo (15 d.p.c., cross-section through the eye) using Goat Anti-Mouse Dkk-1 Antigen Affinity-purified Polyclonal Antibody (Catalog # AF1765) at 15 µg/mL overnight at 4 °C. Tissue was stained using the Anti-Goat HRP-DAB Cell & Tissue Staining Kit (brown; Catalog # CTS008) and counterstained with hematoxylin (blue). View our protocol for [Chromogenic IHC Staining of Frozen Tissue Sections](#).

PREPARATION AND STORAGE

Reconstitution	Reconstitute at 0.2 mg/mL in sterile PBS.
Shipping	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature recommended below. *Small pack size (-SP) is shipped with polar packs. Upon receipt, store it immediately at -20 to -70 °C
Stability & Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. <ul style="list-style-type: none"> ● 12 months from date of receipt, -20 to -70 °C as supplied. ● 1 month, 2 to 8 °C under sterile conditions after reconstitution. ● 6 months, -20 to -70 °C under sterile conditions after reconstitution.

BACKGROUND

Dickkopf related protein 1 (Dkk-1) is the founding member of the Dickkopf family of proteins that includes Dkk-1, -2, -3, -4, and a related protein, Soggy (1, 2). Dkk proteins are secreted proteins that contain two conserved cysteine-rich domains separated by a linker region. Each domain contains ten cysteine residues (1-3). Mature mouse Dkk-1 is a 40 kDa glycosylated protein that shares 86%, 96%, 83% and 82% amino acid (aa) sequence identity with human, rat, rabbit and bovine Dkk-1, respectively. It also shares 41% and 36% aa identity with human Dkk-2 and Dkk-4, respectively. Dkk-1 and Dkk-4 are well documented antagonists of the canonical Wnt signaling pathway (1, 2). This pathway is activated by Wnt engagement of a receptor complex composed of the Frizzled proteins and one of two low-density lipoprotein receptor-related proteins, LRP5 or LRP6 (4). Dkk-1 antagonizes Wnt by forming ternary complexes of LRP5/6 with Kremen1 or Kremen2 (4, 5). Dkk-1/LRP6/Krm2 complex internalization has been shown to downregulate Wnt signaling (4, 5). Dkk-1 is expressed throughout development and antagonizes Wnt-7a during limb development (6, 7). Other sites of expression include developing neurons, hair follicles and the retina of the eye (8, 9). The balance between Wnt signaling and Dkk-1 inhibition is critical for bone formation and homeostasis (10). Insufficient or excess Dkk-1 activity in bone results in increased or decreased bone density, respectively (8, 11). In adults, Dkk-1 is expressed in osteoblasts and osteocytes, and neurons. Cerebral ischemia induces Dkk-1 expression, which contributes to neuronal cell death (12).

References:

1. Glinka, A. *et al.* (1998) *Nature* **391**:357.
2. Niehrs, C. (2006) *Oncogene* **25**:7469.
3. Bullock, C.M. *et al.* (2004) *Mol. Pharmacol.* **65**:582.
4. Mao, B. *et al.* (2001) *Nature* **411**:321.
5. Mao, B. *et al.* (2002) *Nature* **417**:664.
6. Kemp, C. *et al.* (2005) *Dev. Dyn.* **233**:1064.
7. Adamska, M. *et al.* (2004) *Dev. Biol.* **272**:134.
8. Li, J. *et al.* (2006) *Bone* **36**:754.
9. Verani, R. *et al.* (2006) *J. Neurochem.* **101**:242.
10. Pinzone, J.J. *et al.* (2009) *Blood* **113**:517.
11. Morvan, F. *et al.* (2006) *J. Bone Miner. Res.* **21**:934.
12. Cappuccio, I. *et al.* (2005) *J. Neurosci.* **25**:2647.